

Examination & Diagnostic Breast Trainer

Part No: 40044

**Learning Has Never
Been More Lifelike**

Responds to gentle palpation
in an extremely lifelike way

Breast is warm to the touch



Benign cyst
for aspiration

Simulated carcinoma
& fibroadenoma

Lightweight & compact base allows for:
- self-examination
(when held up against chest)
- table top presentation

Skin surface is washable
using soap & water

Variable cyst volume adjusted using
a hidden refillable reservoir attached
to the underside of the breast

Examination & Diagnostic Breast Trainer

Part No: 40044

A life-size realistic breast used for training in examination procedure & diagnostic techniques.

Skills Gained:

- Recognised clinical breast examination technique
- Diagnosis of pathologies:
 - carcinoma
 - fibroadenoma
- Aspiration of a cyst
- Self-examination technique



- Simulated carcinoma & fibroadenoma
- Benign cyst for aspiration

Product Qualities:



- Breast is warm to the touch



- Responds to gentle palpation in an extremely lifelike way



- Variable cyst volume (adjusted using a hidden refillable reservoir attached to the underside of the breast)



- Lightweight & compact base allows for:
 - self-examination (when held up against chest)
 - table top presentation



- Skin surface is washable using soap & water

You'll Receive:

- 40040 1 Examination & Diagnostic Breast
- 40106 1 Breast Aspirant (50ml)
- 1 Rigid Two-Part Plastic Base

Packed weight (kg) & dimensions (mm) 1.25 kg 320 l 270 w 140 h



Learn more at

limbsandthings.com

Limbs & Things Ltd. Sussex Street, St Philips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 F: 00 44 (0)117 311 0501

E: UK Customers - gbsales@limbsandthings.com

E: International Customers - sales@limbsandthings.com